




Rocket[®] Cervix Adaptors for Hysterosalpingography

The Rocket[®] Cervix Adaptor Range provides a reliable, comfortable and cost effective method of performing uterine canalisation for hysterosalpingography and hysterosalpingo-contrast sonography.



- **TWO VERSIONS:** As illustrated and with roller tube clamps to retain contrast media and vacuum.
- **COMPATIBLE** with aqueous and oil based contrast media.
- **LATEX FREE** minimises risk of patient tissue reactions.
- **SMOOTH TIP & LOW PROFILE CUP** ensures low trauma attachment and removal – ideal for the nulliparous cervix.

| Ordering Information | Size | Code |
|--|------------|---------------|
| Rocket[®] Cervix Adaptors sterile, for single patient use. 10 per carton | Small 27mm | R57490 |
| | Large 30mm | R57491 |
| Rocket[®] Cervix Adaptors with tubing clamps sterile, for single patient use. 10 per carton | Small 27mm | R57493 |
| | Large 30mm | R57494 |
| Rocket[®] Cervix Adaptor extension tube sterile for single patient use, packed individually | | R57496 |

