




Rocketmedical

Rocket[®] DUO[™] IUI Catheter



Rocket[®] DUO[™] IUI Catheters are the latest generation of two part insemination catheters designed to improve ease of insertion and minimise uterine trauma.

- **TWO PART CONSTRUCTION IUI** with smooth outer sheath with high memory characteristics ensures ease of insertion, particularly in those patients where access to the external os is difficult.
- **SOFT INNER IUI** with fully domed distal tip and side eyes for optimum insemination
- **COMPETITIVELY PRICED** to provide excellent value in routine IUI.



Specifications:

- (a) Inner Transfer Catheter: Polyurethane OD: 1.4 mm, ID: 0.9 mm, fully domed, dual side eye.
 (b) Outer Sheath: HDPE, 5 x 1cm marking from distal tip.
 (c) Inner Catheter volume: 18cm: 0.5ml

Rocket [®] DUO [™] Range	Length	Code
Rocket[®] DUO[™] IUI Catheters. Sterile, for single patient use, packed in cartons 25 units	18cm	R57625