



Rocketmedical
Medical design and manufacture in the UK

Customer Information Bulletin:

Changes to inner catheter luer connectors on Rocket® Embryo Transfer Catheters

Dear Customer,


As part of our continuous improvement programme, we are introducing a new style moulded luer fitting to the inner catheters of all our embryo transfer and trial transfer catheters.

Traditionally, we have used a 'winged' luer connector on Embryon® Embryo Transfer inner catheters, however, during development of the highly successful **Rocket® BulbTip™** Range of Embryo Transfer Catheters, we introducing a new, easy grip, moulded connector.


From **1st December 2006** we will be rolling out this improved moulded connector throughout the Rocket® ET Range.

The change applies to these versions and custom derivatives:

- Embryon® Embryo Transfer & Trial Transfer Catheters
- Embryon® ThinWall™, ThinWall™ SOFT™ Embryo Transfer & Trial Transfer Catheters
- Rocket® Genesis & Microbore ET Catheters
- Rocket® DUO™ IUI Catheters



Old style



New style, fully moulded

**If you have any further questions, please contact our
Customer Services Team for assistance: 0191 419 6988**

Rocket Medical plc. Sedling Road. Washington NE38 9BZ.
Head Office & Export Customer Services: +44 (0) 192 365 1404

Rev 02: Issued: 03/12/06. Revised 26/04/07