



**Rocketmedical**  
Medical design and manufacture in the UK

# Customer Information Bulletin: Rocket® Indwelling Pleural Catheter Insertion Set

Dear Customer,

As part of our continuous product development program, we are introducing an upgrade to the Rocket® Indwelling Pleural Catheter Insertion Set R54400-16-40.

## New Short Dilator Introduced

16G x 11cm Dilator.



Following concerns by members of the British Thoracic Society, the reduction to 11cm is to minimise the risk of over insertion during catheter placement.

If required the larger 16Gx17cm Split Sheath Dilator is available via order number R54410-16-17 via customer services.

Customer Information Bulletins on all of our products can be found on our web site at:  
[www.rocketmedical.com/customer\\_info\\_main.html](http://www.rocketmedical.com/customer_info_main.html)

We have also introduced a free email based notification service to advise you when new CIB's, relevant to your clinical area are issued. You can subscribe to this free service by going to:  
[www.rocketmedical/cib\\_signup.html](http://www.rocketmedical/cib_signup.html)

**If you have any questions, please contact your nearest Customer Services Team for assistance:**

<b>UK</b>	0191 419 6988
<b>EU, Europe, Middle East, Far East, Australia, NZ.</b>	+44 1923 651404