




Customer Information Bulletin

Rocket® Multipurpose Needle Drain Set




R58800-08-ND Rocket® Needle drain with new straightener sleeve.

Dear Customer,

As part of our continuous product improvement programme, we are introducing a minor change to the R58800-08-ND Rocket® Multipurpose Needle Drain.

The Rocket® Multipurpose Needle Drain Set offers a purpose designed, cost-effective solution for routine fluid drainage procedures and is specifically designed for the drainage of abdominal ascites. A key feature of the design is the protective pigtail tip which encloses the drainage holes, minimising the risk of occlusion from contact omentum or other abdominal structures.

The Rocket® Needle Drain has an EasyTip™ core which prevents the soft catheter material from ruckling back during placement, so improving ease of insertion and minimising patient discomfort.



Change:

To improve ease of loading the needle, a new design red coloured outer sleeve has been introduced.

Simply slide the sleeve forward to fully unravel the pigtail and slide the needle completely into the drain. The luer fittings will then securely lock together with a half anticlockwise turn.



Always remove the red sleeve before insertion

Customer Information Bulletins on all of our products can be found on our web site at:

www.rocketmedical.com/customer_info_main.html

We have also introduced a free email based notification service to advise you when new CIB's, relevant to your clinical area are issued. You can subscribe to this free service by going to:

www.rocketmedical/cib_signup.html

If you have any questions, please contact your nearest Customer Services Team for assistance:

**In the UK - call Customer Services:
0191 419 6988**

**Outside of UK - call Customer Services:
+44 192 365 1404**