



# Customer Information Bulletin:

## Rocket® Fetal Blood Sampling Kits improved contents, added value.

**Rocket® FBS Kits** have previously contained either 6 x cotton swab sticks or 5 x gauze swabs and a metal 'Rampley' style sponger holder. This latter version was produced for those few customers who preferred the greater rigidity of the metal sponge-holder.

In response to customer demand, we have looked at a number of ways of combining these two versions to give all our customers the benefits of both types of kit. Now our colleagues at Pelican Healthcare have produced an excellent, single use plastic sponge-holder which has the required rigidity and quality to meet our needs and those of our customers.

As a result, with effect from **1<sup>st</sup> February 2007** we will be adding 5 x gauze swabs and a single use sponger holder into every Rocket® FBS Kit.

The two kits, **R57012-R** and **R57013-R** containing the old style metal 'Rampley' sponge-holder will now be discontinued.

- **NEW AMNILUME™ COMPATIBLE:** all Rocket® FBS Kits are designed for use with the revolutionary Amnilume™ Single Use Light Source system.
- **FOUR COMPREHENSIVE KITS:** each designed to complete up to 5 individual FBS samplings over the delivery period. Each kit contains:

- |   |                                    |
|---|------------------------------------|
| ✓ Endoscope & Obturator compatible with Amnilume™ Single Use Light Source | ✓ 1 x Rocket® FBS Blade Holder     |
| ✓ 6 x 12" cotton swab sticks  | ✓ 1 x Yellow Soft Paraffin BP      |
| ✓ 5 x gauze swabs   | ✓ 1 x capillary tube holder        |
| ✓ 1 x Disposable sponge holder  | ✓ 10 x capillary tube stoppers     |
| ✓ 2 x Rocket® FBS Blades  | ✓ Tray wrapped in waterproof drape |

Packed individually in cartons of 10 units

Additional contents	R57012	R57013	R57018	R57024
3 x long heparinised tubes, glass, straight	✓			
3 x long heparinised tubes, glass, angled		✓		
5 x 90µl plastic Rocket® capillary tubes	✓	✓		
5 x 55µl plastic RADIOMETER® safeCLINITUBES			✓	
5 x 85µl plastic RADIOMETER® safeCLINITUBES				✓
Compatible with <b>NEW</b> Amnilume™ Single Use Light Source	✓	✓	✓	✓

**If you have any further questions, please contact your nearest Customer Services Team for assistance:**

<b>UK</b>	<b>0191 419 6988</b>
<b>EU, Europe, Middle East, Far East, Australia, NZ.</b>	<b>+44 1923 651404</b>