

# Rocket® Ascitex™

## Multipurpose Drainage Procedure Pack

**Rocket Medical plc** have developed a cost-effective solution to the tiresome problem of assembling a wide range of ward consumables, often at short notice, to complete a routine drainage procedure.

The **Rocket® Ascitex™ Multipurpose Drainage Procedure Tray** contains all the items necessary to safely and conveniently complete the procedure in the department or ward, using the long established Seldinger insertion technique.



- **MULTIPURPOSE DRAIN** with soft, kink resistant catheter material, improves patient comfort and prevents unnecessary catheter replacement caused by harder materials.
- **NITINOL GUIDEWIRE** with soft, non traumatic platinum tip to ensure safe delivery of the catheter into the abdominal cavity or organ.
- **CHOICE** of classic Seldinger or marked Chiba insertion needles
- **SECURE LOCKING SYSTEM** to give reliable pigtail locking without cumbersome 'keys' or bulky hubs. Self adhesive removal instructions (for addition to the patient's notes) are provided with every locking drain.
- **COMPLETE PROCEDURE TRAY** including syringes, needles, drapes, swabs and ascites drainage bag with tap

**Suitable for:**

- Abdominal paracentesis
- Pneumothorax, spontaneous or iatrogenic.
- Effusions, benign or malignant, except those with high protein content.
- Low viscosity empyema.
- Nephrostomy

Description	Code
<p><b>Rocket® Ascitex™ Multipurpose Drainage Procedure Pack</b></p> <p>8.4FG locking pigtail catheter, 6 holes x 30cm, 20ml, 12ml syringes, leur lock, 21g x 1.5" needle, 25g x 5/8" needle, Scalpel No.11, 18G x 15cm Chiba Needle, 18G x 7cm Seldinger Needle, Nitinol guidewire 0.035" x 80cm x 6cm soft tip, 9fr x 8in dilator, 30cm drain bag connection set + 2 way tap, 2000ml drainage bag with drain port., 5in Sharp/Blunt Scissors, 60ml gallipot, <i>Disposafe</i> sharps mat, Gauze swabs x 5, Sponge stick, 8", 11x15cm <i>Mepore</i> dressing, 50x60m adhesive fenestrated drape, Insertion and removal instructions, 100 x 100cm sterile field wrap. Sterile, for single patient use, supplied individually</p>	<p><b>R58800-84-30</b></p>

