





Rocket® Coagulator Probe Cover

The Rocket® Coagulator Probe Covers offer a simple, cost effective solution for the prevention of microbial contamination related to the reprocessing of autoclavable and non autoclavable WISAP® Coagulator Probes.

- **OPERATING TECHNIQUE IS UNAFFECTED.** Users should continue to use their normal treatment protocols and repeat if required to achieve a satisfactory post treatment appearance. Rocket Medical normally recommend probe application of 100°C for 30 secs.
- **EASY** 'roll-on' application with pre-lubricated sheath. Simple removal and disposal for busy clinic usage.
- **BIOCOMPATIBLE USP Class VI silicone** guaranteed 100% latex free ensures minimal patient tissue reaction.
- **DESIGNED & TESTED** to BS EN ISO 4074:2002 to ensure a reliable microbial barrier to eliminate cross infection risks from coagulator probes
- **STERILE, SINGLE USE** eliminates costly reprocessing concerns
- **COST EFFECTIVE** to reduce patient costs
- **COMPATIBLE** with new style WISAP® coagulator probes types: A,B,C,F,G & H.

Probe Choices:



Description	Order Code
Rocket® Coagulator Probe Covers for use with WISAP® Coagulator probes types: A,B,C,F,G & H. For single patient use, sterile, supplied in cartons of 20 units.	R57470

WISAP® is a registered trademark of WISAP GmbH.

