


Rocket® Bone Marrow Aspiration Needles

The **NEW Rocket®** Bone Marrow Needle is the direct result of recent user research. Our new design includes those features which clinicians most commonly need or desire in the ideal device for bone marrow aspiration.

- **SECURE 'T' GRIP HANDLE** which sits comfortably in the operator's hand for confident control of the needle during insertion.
- **DEPTH CONTROL SLEEVE** gives total control of depth of penetration allowing the user to set needle lengths between 25-60mm for safe sternal and iliac crest puncture.
- **15G NEEDLE** with precision ground tip to enable easy penetration of bone.



Ordering Information	Code
<p>Rocket® Bone Marrow Aspiration Needles 15g needle with ergonomic moulded 'T' handle, retaining plate, stylet and removable sleeve variable between 25-60mm needle protrusion. Peel pouch packed in cartons of 10 needle sets, sterile, for single patient use</p>	<p>R56780-00-T</p>

