

Rocket®

Amniocentesis & CVS Sampling Set



The Rocket® Amniocentesis & CVS Set has been developed with the assistance of a leading UK Fetal Medicine Unit and is the result of a program of continuous product development and design improvement.

- **HIGH RESOLUTION ECHO BEVEL DESIGN** to give smoother, easier tissue penetration with pinpoint ultrasound definition of the needle tip in either amniocentesis and CVS.
- **NEEDLE MARKINGS** and depth marker help ensure correct positioning and minimise the risk of over insertion.
- **INTEGRAL TUBING SET** – purpose designed 10cm extension tube set is included in every pack for simple, safe attachment to the needle without the need to find costly off-the-shelf items which cannot offer guaranteed compatibility and security.
- **SECURE LUER LOCK FITTINGS** ensure the needle, extension set and syringes remain correctly and securely connected throughout
- **CE MARKED, PURPOSE DESIGNED DEVICE**, designed specifically for amniocentesis and transabdominal chorionic villus sampling which avoids the risks of inappropriate use of non specific generic devices.



Description	Size	Code
Rocket® Amniocentesis & CVS Sampling Set echogenic amniocentesis needle with 1cm markings and 10cm luer lock extension tube set. Sterile for single use, packed in cartons of 10 sets	20G x 15cm	R57036

