

Rocket®

Amniocentesis Exchange Set



The Rocket® Amniocentesis Exchange Set has been developed with the assistance of a leading UK Fetal Medicine Unit and is the result of a program of continuous product development and design improvement.

- **HIGH RESOLUTION ECHO BEVEL DESIGN** to give smoother, easier tissue penetration with pinpoint ultrasound definition of the needle tip in either amniocentesis and CVS.
- **NEEDLE MARKINGS** and depth marker help ensure correct positioning and minimise the risk of over insertion.
- **INTEGRAL TUBING SET** – purpose designed 10cm extension tube set is included in every pack for simple, safe attachment to the needle and supplied syringe without the need to find costly off-the-shelf items which cannot offer guaranteed compatibility and security.
- **10ml 3 RING SYRINGE & 3 WAY STOPCOCK INCLUDED:** - each set contains a 10ml, 3 ring syringe and widebore 3 way stopcock for comfort and security when repeated aspiration/infusion is required.
- **SECURE LUER LOCK FITTINGS** ensure the needle, extension set stopcock and syringe remain correctly and securely connected throughout the procedure.



Description	Size	Code
Rocket® Amniocentesis Exchange Set echogenic amniocentesis needle with 1cm markings, 10cm luer lock extension tube set, 3 way- stopcock, 3 ring syringe. Sterile for single use, packed in cartons of 10 sets	20G x 15cm	R57035

