



# Rocket®

## Amniotomy Hooks

Rocket Medical plc have a long and enviable reputation for reliability and innovation in the Labour Ward and our acquisition of the Copeland® range complements our existing fetal and labour ward products. As a result, we can now offer two versions to meet your individual needs

**Both Rocket® Amniotomy Hooks feature:**

- **ROUNDED DISTAL END** to prevent damage during insertion
- **GUARDED HOOK:** located at a controlled depth to permit perforation of the membranes without trauma to the surrounding tissue. (See image detail)
- **STRAIGHT: SURGI-HOOK™** with contoured handle to indicate angle of perforating point, Hook faces away from direction of insertion to protect against accidental tear.
- **CURVED:** anatomically shaped to ensure a positive grip in the users hand and to follow the normal curved pelvic outlet. for precise finger location and positioning.



Description	Per Carton	Code
<b>Rocket® Amniotomy Hook, curved</b> 25cm. Sterile, for single patient use, individually packed,	<b>100</b>	<b>R57065</b>
<b>Rocket® Surgi-Hook, straight</b> 26cm. Sterile, for single patient use, individually packed,	<b>100</b>	<b>R57066</b>

