






Rocket FBS Accessories

The Rocket FBS Accessories are purpose designed to complement the Rocket FBS Kit Range providing cost effective flexibility to meet your patient's needs.

Description		Code
Rocket Amnilume™ Disposable Light Source compatible with all Rocket Single Use FBS Amnioscopes Cartons of 10 units, sterile, for single patient use		R57004
Rocket White Soft Paraffin BP (Petroleum Jelly) 3g tube. Carton 20 units, sterile for single patient use		R57006
Rocket FBS Blades. 2.2mm, angled. Cartons 100 units, sterile for single patient use		Plug: R57010 Screw: R57011
Rocket FBS Blade Holder & Blade. 2.2mm screw type, angled. Cartons 10 units, sterile for single patient use		R57011-HLDR
Rocket Long Glass Heparinised Tubes 295 x 3.0 x 1.0mm 100 IU/ml Sodium Heparin. Carton 100 tubes, non sterile	Straight Angled	R57014 R57015
Rocket Amnioscope complete with obturator Compatible with Amnilume™ Disposable Light Source (not included) 142mm long x 33mm distal OD, 52mm proximal OD 200mm long x 33mm distal OD, 52mm proximal OD Cartons of 10 sets, sterile, for single patient use		R57016 R57016-PLUS
Rocket Long Amnioscope without obturator, complete with Amnilume™ Disposable Light Source 200mm long x 23mm distal OD, 52mm proximal OD Cartons of 10 sets, sterile, for single patient use		R57017
Rocket Precision Capillary Tubes – Plastic, High Heparin Carton of 10 packs x 5 capillary tubes, sterile for single patient use	55µl 85µl 100µl	R57018-TUBE R57024-TUBE R57026-TUBE

