


# Rocket

## Disposable FBS Endoscopes



Rocket Medical is unique in producing a range of THREE Amnioscopes/Endoscopes to meet the widest range of patient needs. All designed to be used with the revolutionary Amnilume™ light source to provide clear, un-obstructed, brightly lit access when time is precious.

- **THREE SIZES** to meet widest range of patient needs:
  - **Standard:** 142mm long x 33mm distal OD with obturator
  - **Early:** 200mm long x 23mm distal OD without obturator, includes Amnilume™
  - **Plus:** 200mm long x 33mm distal OD with obturator
- **COMPATIBLE** with Rocket Amnilume™ disposable Light Source
- **COVENIENT** and **COST EFFECTIVE:** Packed sterile, for single patient use in cartons of 10 units, means that you can easily accommodate your patient's individual needs without storing bulky and costly stock when space and budgets are limited.



Ordering Information		Code
<b>Rocket FBS Endoscope (Amnioscope)</b> in cartons of 10 units, sterile, for single patient use <b>Standard:</b> 142mm long x 33mm distal OD with obturator <b>Plus:</b> 200mm long x 33mm distal OD with obturator	<b>Standard</b>	<b>R57016</b>
	<b>Plus</b>	<b>R57016-PLUS</b>
<b>Rocket FBS Endoscope (Amnioscope) complete with Amnilume Light Source</b> in cartons of 10 units, sterile, for single patient use <b>Early:</b> 142mm long x 23mm distal OD without obturator	<b>Early</b>	<b>R57017</b>

