

Rocket®

Fetal Blood Sampling Kits

In conjunction with the world's leading manufacturers of blood gas analysers, Rocket Medical have specifically developed kits offering reliable, safe and cost effective solutions for routine fetal blood sampling.

■ **RELIABLE, COMPATIBLE ANALYSIS:** Rocket® FBS Kits are specifically designed for use with these systems:

- Siemens® (Bayer) RapidLab™ Series
- Radiometer® ABL™ 70, 90, 800, 700, 77, and 5/500 Series
- Instrumentation Laboratory® GEM™ 3000 & 4000 systems
- ROCHE micro volume systems.

■ **COMPREHENSIVE KITS:** each designed to complete up to 5 individual FBS samplings over 30 mins. Each kit is ready for immediate use to minimise delay when time is precious.

■ **MAXIMUM ILLUMINATION** with the minimum fuss, all Rocket® FBS Kits now come with the revolutionary Rocket® Amnilume™ Single Use Light Source which produces 8800mcd @ 30°.

■ **X-RAY DETECTABLE SWABS:** To meet NPSA recommendations following publication of NPSA Rapid Response Report NPSA/2010/RRR012: 'Reducing the risk of retained swabs after vaginal birth and perineal suturing'. all Rocket FBS Kits now contain green, x-ray detectable swabs.

- | | |
|-------------------------------------|------------------------------------|
| ✓ Endoscope & Obturator | ✓ 1 x Rocket® FBS Blade Holder |
| ✓ Amnilume™ Single Use Light Source | ✓ 1 x Soft Paraffin BP |
| ✓ 6 x 12" cotton swab sticks | ✓ 1 x capillary tube holder |
| ✓ 5 x gauze swabs | ✓ 10 x capillary tube stoppers |
| ✓ 1 x disposable sponge holder | ✓ Tray wrapped in waterproof drape |
| ✓ 2 x Rocket® FBS Blades | |



Compatible Analyser	Recommended sample volume	Rocket® FBS Kit inc Amnilume™	5 x extra sample tubes
ABL™ 5/500, 70, 77 & 90 Series	>65µl	R57024	R57024-TUBE
ABL™ 800 FLEX with FLEXMODE selected	>35µl	R57012 R57013 R57018 R57024	R57018-TUBE R57024-TUBE
ABL™ 700 Series with tube volume selected	>35µl	R57018	R57018-TUBE
Instrumentation Laboratory (IL) GEM™ 3000 Series	170µl	R57027	
Instrumentation Laboratory (IL) GEM™ 4000 Series	65µl	R57024	R57024-TUBE
ROCHE® Systems	55µl	R57018	R57018-TUBE
SIEMENS® (Bayer) Rapidlab™ Series	100µl	R57026	R57026-TUBE

