

# Rocket® Pudendal & Paracervical Block Needle



Rocket® Pudendal Block Needle with red safety clip – remove before use

**Rocket® Pudendal Needles provide a safe and controllable method for delivering a Pudendal nerve block.** This route typically provides satisfactory perineal anaesthesia for normal delivery, low forceps manipulation and episiotomy. The advantage of this approach over a major conduction block is the lack of circulatory and respiratory disturbance. Typically the course and duration of labour are unaffected, even though there may be some impairment of the perineal-diaphragmatic-abdominal reflex.

- **LUER LOCK NEEDLE** provides up to a maximum of 10mm tissue penetration for increased safety.
- **REMOVABLE FINGER RING** can be re-sited at Position 1 or left as supplied in Position 2 to give best individual control and comfort during use.



Description	Code
<b>Rocket® Pudendal &amp; Paracervical Block Needle 20G x 130mm</b> Sterile, for single patient use, individually packed, 20 per carton	<b>R57040</b>

