

Rocket®

Chest Drainage Solutions

The Rocket® Blue Bottle™ range is the UK's biggest selling, underwater seal chest drainage system.

For a complete listing of all our chest drains and tubing set combinations, please visit our website: www.rocketmedical.com or contact our Customer Services Team to arrange a demonstration from one of our Clinical Sales Specialists.

The system features:

- **UNDERWATER SEAL CHEST DRAINAGE BOTTLE**
for either single collection (R54500) or double collection (R54509)
- **UPRIGHT, INHERENTLY STABLE DESIGN**
which does not require tiresome stands or carriers
- **1800ml CAPACITY:**
available in two graduations 25ml: R54500 or 5ml: R54509
- **AVAILABLE WITH A WIDE RANGE OF TUBING SETS:**
including, single cone connector for connection to large bore catheters: (R54502) and also for use with Rocket® Seldinger Chest Drainage Sets: (R54539)
- **SIMPLE ORDERING:**
all products are available direct, many through NHS Logistics.
NHS Codes listed below.



Description	Code	NHS Code
Rocket® Blue Bottle™ Single chamber, 1800ml capacity, sterile, 10 units per carton.	R54.500	FSW076
Rocket® Double Chamber: ideal for paediatric or specific drainage recording. 1800ml capacity bottle, 1 st chamber: 5ml, 2 nd chamber: 10ml graduations, sterile, 10 units per carton	R54.509	FSW047
Rocket® Tube Set: with bottle cap, integral straw, support spring, 1.5mm with fir tree connector, sterile, single wrapped, 10 units per carton	R54502	FSW077
Rocket® Tube Set with bottle cap, integral straw, support spring, 1.5mm with plain end, suitable for connection to Rocket® Seldinger Chest Drain Sets, sterile, single wrapped, 10 units per carton	R54539	FSW187

