

# Rocket®

## Pleuraseal® III – Triple Bottle System

The Rocket® Pleuraseal® III water seal chest drainage unit is a three chamber system comprising a suction control chamber, water seal chamber and a fluid collection chamber.

For further information please visit our website: [www.rocketmedical.com](http://www.rocketmedical.com) or contact our Customer Services Team to arrange a demonstration from one of our Clinical Sales Specialists.

The system is indicated for:

- Areas where accurate suction control and drainage measurement are required.

The system features:

- Easy to prime water seal chamber via attached funnel.
- Height of the water in the suction control chamber allows for the prescription of desired suction levels.
- Fluid collection chamber with 2400 ml capacity.
- Positive pressure release valve.
- High negativity float valve and filtered manual high negativity vent.
- All products are available direct from our Customer Services Team on **0191 419 4488**



Description	Code	NHS Logistics Code
<b>Rocket® Pleuraseal® III Dual Collection System</b> Dual chamber, 2400ml capacity, 1ml-5ml graduations, 1.5m patient tube set, sterile, 5 units per carton.	<b>R54542</b>	
<b>Rocket® Pleuraseal® III Single Collection System</b> Single chamber, 2400ml capacity, 1ml-5ml graduations, 1.5m patient tube set, sterile, 5 units per carton.	<b>R54543</b>	<b>FSW371</b>
<b>Rocket® Pleuraseal® III Single Collection System with short patient connection tube</b> Single chamber, 2400ml capacity, 1ml-5ml graduations, 1m patient tube set, sterile, 5 units per carton	<b>R54545</b>	

