

Rocket® Craft™ Pump Annual Servicing Package

Authorised, regular servicing of your Rocket® Craft™ Pump is vital to maintain its long term performance and reliability. To make servicing even easier, our new Annual Servicing Package provides a fixed price for the standard annual maintenance which each pump requires.

Rocket® Craft™ Pumps have a long and enviable reputation for reliability and typically only require a standard service programme and parts replacement to keep them in first class order. Where a pump requires additional components such as a replacement casing following accidental damage, these will be repaired at the time of service at an additional charge.

- **FACTORY SERVICING:** Where better to have your Rocket® Craft™ Pump serviced than in the same factory and by the same technicians that first built it?
- **FIXED PRICE ANNUAL SERVICE PACKAGE:** Includes factory inspection, replacement of service pack items, lubrication, adjustment and re-calibration.
- **PARTS INCLUDED:** These items are replaced as standard and are included in the service cost:
 - Intake filters x 2
 - Pump diaphragm
 - Diaphragm top plate
 - Valve set
 - Gasket set
- **NOT INCLUDED:** Parts and labour to replace items not covered by the Annual Service pack such as damaged casings and broken gauges. Items incurring an additional charge will be notified to the customer before work is completed. Return delivery charges are not included in the Annual Servicing Package.
- **TO ARRANGE A SERVICE:** Contact our Customer Services Team on +44 (0) 191 419 6989 or email us at customerservices@rocketmedical.com



Ordering Information

Rocket® Annual Servicing Package
suitable for all Craft™ Suction Pumps. Delivery charges apply.

R29670-SERV

