


Rocket® DUO™ IUI Catheter



Rocket® DUO™ IUI Catheters are the latest generation of two part insemination catheters designed to improve ease of insertion and minimise uterine trauma.

- **NEW: TWO PART CONSTRUCTION IUI** with smooth outer sheath with high memory characteristics ensures ease of insertion, particularly in those patients where access to the external os is difficult.
- **NEW: SOFT INNER IUI** with fully domed distal tip and side eyes for optimum insemination
- **COMPETITIVELY PRICED** to provide excellent value in routine IUI.
- **STERILE** - Packed in peel pouches, 25 per carton, intended for one time use.



Specifications:

- (a) Inner Transfer Catheter: Polyurethane OD: 1.4 mm, ID: 0.9 mm, fully domed, dual side eye.
- (b) Outer Sheath: HDPE, 5 x 1cm marking from distal tip.
- (c) Inner Catheter volume: 18cm: 0.5ml

Rocket® DUO™ Range	Length	Code
Rocket® DUO™ IUI Catheters. Sterile, for single patient use, packed in cartons 25 units	18cm	R57625

